

OURS GOES TO 10.



PerkinElmer ExacTag

Enabling the highest multiplexed quantitative protein expression profiling

PerkinElmer ExacTag™ kits for protein expression profiling enable simultaneous quantification of up to 10 samples using any tandem mass spectrometer. Label samples at the protein level and then pool for fractionation and digestion to allow for:

- sample complexity reduction
 - low abundance protein analysis
 - improved selectivity
- ...all without sacrificing precision.

- 2, 4, 7 and 10 multiplex levels available, for labeling either cysteine or lysine residues
- High level multiplexing allows accurate pharmacological profiling from protein expression data
- Protein labeling early in workflow improves quantitative precision
- Preparation and analysis performed on pooled samples for fast, efficient workflow
- Compatible with any tandem mass spectrometer
- Compatible with all standard methods of sample complexity reduction including gel electrophoresis, LC & immunoprecipitation



Visit our web site for a special introductory offer on ExacTag

www.perkinelmer.com/gototen



www.perkinelmer.com/gototen